



P090VL80

P090VL81

STERILE R

PETRI DISH



P090VL80

P090VL81

STERILE R



IVD

STERILE R



Material:

Polistirene

Model:

90 x 14 mm

Standard packaging:

Bags of **24** pcs

Box of **600** pcs

∅	External dish:	89,14 mm
∅	Internal dish:	86,40 mm
	Dish height:	13,80 mm
∅	External lid:	92,90 mm
∅	Internal lid:	90,70 mm
	Lid height:	7,00 mm
	Assembled height:	14,80 mm
	Total weight:	13,70 g



Ventilation_

The ventilated plates have **8** vents, **4** of which are functional for the passage of air and **4** additional to guarantee maximum stability of the lid.



Stability_

Some improvements have been adopted to guarantee maximum stability of the dishes once stacked, including:

- **Increased** plastic thickness for greater resistance to deformation under thermal and mechanical stress
- **4 vents** that help center the plate and lid
- **Oversized rings** on the edge of the lid and on the bottom of the plate to further improve its stability.

Possibility of ordering the sterile product in a triple bag, in compliance with pharmaceutical production standards, ref. P090VL81TS.



anicrin srl

Via Moglianese, 44/L
30037 Scorzé (VE) Italy

Tel.: + 39 041 445 660

info@anicrin.com
www.anicrin.com

A brand of:
i'M